

We've been hearing good things...



"I was so amazed at the instant results that I had to have one immediately. I am extremely grateful to you for this product that has increased my quality of life quite a bit."

— Eric G., Eden Prairie MN

"Used my EarPopper before scuba diving — no problems descending. No pressure build-up in my head — a much smoother dive."

— Kerri B., Scottsville KY

"I can now relax and enjoy flights without worrying about extreme ear pain when descending for landing."

— Barbara R., Solena Beach CA

"After searching for other options online, we found your website. We asked our pediatrician to prescribe the EarPopper. Our son just went back for a check-up 6 weeks later... COMPLETELY FLUID FREE!"

— Abby N., Bloomington IN

"Within three weeks of use my daughter's hearing improved remarkably."

— J Kevin B., Moore OK

"Helped unclog my Eustachian Tubes better than any other method and without meds!!"

— Claudia B., Delano MN

How to get the EarPopper:

1. Obtain a prescription from your healthcare provider for an EarPopper.
2. Take your prescription along with this brochure to any pharmacy.
3. The pharmacy will order the EarPopper and should have it available within 24 hours.

PHARMACISTS: *The EarPopper can be ordered through your medical supply company.*

National Drug Code:	08564 000000
McKesson Item Number:	1777259 EARPOPPER
AmerisourceBergen Item Number:	762-551 EARPOPPER
Cardinal Item Number:	4081691 EARPOPPER

If you have any problems, please contact Micromedics Customer Service at 800-624-5662 or visit www.earpopper.com



In the U.S.A., the EarPopper is FDA cleared and is available by prescription only.



EarPopper is CE registered in the EU.



PHONE: 651 452-1977 | FAX: 651 452-1787
1270 EAGAN INDUSTRIAL RD | ST. PAUL, MN 55121

www.earpopper.com



**No Drugs!
No Ear Tube
Surgery!**

**Ear Pressure
Relief Device**



The EarPopper

For ear pressure relief!

Everyone can relate to the discomfort of ear pressure experienced while flying or riding a fast moving elevator — especially when combined with a cold or allergies. Pressure on the eardrums can be painful, cause hearing loss and lead to long-term problems if the pressure isn't relieved.



The EarPopper treats ear problems caused by:

- Flying
- Scuba Diving
- Allergies
- Colds
- And much more

The EarPopper treats the problem at the source — The Eustachian Tube

The Eustachian Tube (ET) is the body's only mechanism for equalizing pressure and ventilating the space behind the eardrum (Middle Ear). The ET normally opens when you swallow or yawn. A normally-functioning ET keeps the middle ear ventilated and the pressure equalized.

When the ET doesn't function normally, the middle ear becomes a closed, unbalanced space which forces uncomfortable pressure on the eardrum while affecting your hearing. This closed space cannot drain properly — resulting in fluid build up behind the eardrum.

The EarPopper balances pressure in the middle ear by delivering a safe, constant stream of air into the nasal cavity. During the moment of swallowing, the air opens the Eustachian Tube, clearing the middle ear, relieving negative ear pressure and allowing any fluids to drain.

EarPopper — the alternative to medication and ear tube surgery

Ask your doctor about a prescription today!



OPEN EUSTACHIAN TUBE



BLOCKED EUSTACHIAN TUBE



Using The EarPopper: As quick as a swallow!

1

Hold nosepiece firmly against nostril opening — creating a tight seal. Plug the other nostril closed.

2

Push button to start the airflow and swallow while the device is running.

3

Repeat on other nostril. After 5 minutes, repeat steps 1-3. This will complete one treatment.

EarPopper